

## FUM-750 Peninsular Fume Hood for Group Teaching

Superstructure has sliding sash on two faces and viewing window at end. Posts and sash with safety glass are of stainless steel. Exterior is of electro zinc plated steel with an acid resistant finish. The illuminated interior has adjustable baffles. The raised edge stainless steel top has integral cylindrical sink. Accessories included blower switch with pilot light, a light switch, two 110 volt AC outlets and a cold water and gas combination fixture. Base cabinet is sold separately. Overall size including base cabinet is 59" long, 30" wide, and 7' – 8" high.

