

## PERCHLORIC ACID FUME HOODS

This superstructure is designed to minimize the potential hazards of activities involved and organic materials are not used in its interior. The entire interior, with easily cleaned covered corners, is of Type 316 stainless steel. A remote controlled spray above and behind the baffles is provided for wash down procedures, the water falling into the trough at the rear of the working top. The posts and sash with safety glass are also of stainless steel, and the exterior is of electro zinc plated steel with an acid resistant finish. Accessories include two 110 volt AC outlets, a blower switch, a light switch for vapor proof light, and remote control valves for cold water, gas and wash down water.

Base cabinet is sold separately. Overall Dimensions including base understructure are 47" long, 34" wide, and 7'-2" high

